Chromaclone® PVS Impression Material

Description/Indications:

Chromaclone® PVS Impression Material is a superhydrophilic vinyl polysiloxane impression material and is an exciting new advancement in prosthodontics. Chromaclone PVS can be utilized for all final impressions for fixed or removable prosthodontics including crowns, bridges, dentures, partial dentures, and implants or whenever impression accuracy is of the prime concern. Chromaclone PVS putty or Clone Bite can be used for all bite registrations.

Procedure:

CARTRIDGES:

1. Load the cartridge in the dispensing gun.
2. Immediately prior to taking the impression remove the cartridge cap and discard.
3. Gently squeeze the dispenser to bleed enough material from the cartridge to ensure even flow of the base and catalyst. (About 1/4 inch of material).
4. Attach the mixing tip.

5. Squeeze the dispenser in order to extrude the product.
6. After use, store cartridges on end with the used mixing tips attached to protect the remaining materials from contamination.
7. Always bleed the cartridge before each use (see number 3 above).

Use Ultradrant's Impression Mixing Tips attached to impression cartridge with intraoral tip to mix material and inject around preparation or use Ultradrant's Impression FX Tip and impression syringe to inject and burnish material around preparation.

TIPS:
Use Ultradrant's Impression Mixing Tips attached to impression cartridge with intraoral tip to mix material and inject around preparation or use Ultradrant's Impression FX Tip and impression syringe to inject and burnish material around preparation.

PLEASE NOTE:

The temperature of Chromaclone PVS Impression Material affects the working and setting times.

Ideally it should be kept at room temperature (18-25° C) or (65-77° F).

Avoid contaminating Chromaclone PVS Impression Material with chemicals that could inhibit the set (e.g. acrylic and methacrylate residues, latex rubber and sulphur compounds).

ULTRADRANT UNIVERSAL TRAY ADHESIVE:

Use of the tray adhesive is very important. For an accurate impression, coat tray with a thin layer including the edge.

Let adhesive dry for 2 minutes before making impression.

TRAYS:

We recommend Chromaclone plastic disposable trays or any other rigid trays. We recommend applying a thin layer of tray adhesive to all trays.

DISPENSING AND GUN ASSEMBLY:

The Ultradrant Impression Gun can be sterilized using heat sterilizers, including STATIM Autoclaves, conventional autoclaves, chemical vapor, or dry heat sterilizers.

DISINFECTION:

The impressions should be rinsed under cold running water after removal from the mouth. Soaking for 30 minutes in liquid disinfectants such as 2% glutaraldehyde (acidic), 0.5% - 1% sodium hypochlorite, 0.1% - 1.5% providone-iodone (halogenated phenol) or 0.13% glutaraldehyde (neutral) will not affect the surface or dimensions.

CASTING:

The impression may be poured immediately after removal from the mouth following disinfection. Surface tension will be reduced if the impression is briefly washed out with a surfactant. Rinse thoroughly afterwards with water.

Ultradrant Impression Materials are compatible with all dental stones, epoxy resins, and polyurethanes.

CONTRA-INDICATIONS

There are no known contra-indications provided when the product is used as directed. There are no known harmful reactions or secondary effects on patients and/or dental personnel using vinyl polysiloxane impression materials.

IMPORTANT NOTES AND PRECAUTIONS

The temperature of Chromaclone PVS Impression Material affects working time. Ideally it should be kept at room temperature between 5-27° C (41-80° F). Even within those boundaries, changes in temperature will affect the working time: A 2° C temperature increase will shorten working time 10-15 seconds; conversely, a 2° C temperature decrease will lengthen working time 10-15 seconds. You may refrigerate to increase working and setting times. Always "bleed" the cartridge before each use by extruding a quarter inch of material before applying a fresh mixing tip. This will ensure even dispensing of fresh base and catalyst pastes. For even more precise application, Ultradrant's Intraoral tips are available.

1. Certain chemical agents like superoxol or hydrogen peroxide may interfere with the set of Chromaclone PVS, care should be taken to remove these residual materials from impression site. Latex gloves (or anything that contains sulfur) and hand lotions can inhibit or contaminate setting of all Chromaclone PVS impression materials.

Powder free or polyethylene gloves are recommended or wash hands thoroughly and mix putty bare handed.

2. Proper care should always be taken to ensure good preparation, cleanliness, and moisture control including hemorrhage and saliva in order to achieve an adequate impression.

3. There is no bubbling or gasing with Chromaclone PVS.

4. Disinfection can be accomplished with any appropriate disinfectant following the original manufacturer's procedures.

5. Multiple casts can be poured from an impression and Chromaclone PVS is compatible with all gypsum products, electroplating and refractory casts for up to 15 days.

6. Temperature affects all set times. All testing for working and setting times was done at 72° F/22° C.

Shelf life is at least 3 years.


7. Do not wear latex gloves or hand lotions while mixing the putty. Wash hands and mix equal amounts of catalyst and base.

COMPATIBILITY

Chromaclone PVS Impression Material is addition-reaction type silicone impression material. It may not be combined with standard or condensation type silicone materials or poly ether materials.

Technical Data:

Heavy, Medium, Light/Wash, Super Light Regular Set	Working Time	2:30 min
	Total Setting Time	5:00 min
Heavy, Medium, Light/Wash, Super Light Fast Set Total	Working Time	1:15 min
	Total Setting Time	2:45 min

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

HAZARD RATING
 4 = Severe
 3 = Serious
 2 = Moderate
 1 = Slight
 0 = Minimal

Chromaclone® PVS
 Impression Material

CE

Manufactured for:
 Ultradent Products Inc
 505 West Ultradent Drive (10200 South)
 South Jordan, Utah 84095 USA
 Made in Italy

EC REP
 Ultradent Products GmbH
 Am Westshover Berg 30
 51149 Cologne Germany

ULTRADENT PRODUCTS, INC.

For product SDS please see our website: www.ultradent.com
 1-800-552-6512; 801-572-4200

EN- See Instructions
 DE- Siehe Anleitung
 FR- Cf. instructions
 NL- Zie gebruiksaanwijzing
 IT- Vedere le istruzioni
 ES- Vea las instrucciones
 PT- Veja as instruções
 SV- Se anvisningar
 DA- Brug anvisningen
 FI- Lue käyttöohjeet
 EL- Άναγνώστε τις οδηγίες

EN- Lot Number
 DE- Chargennummer
 FR- Numéro du lot
 NL- Batchnummer
 IT- Numero del lotto
 ES- Número de lote
 PT- Número do lote
 SV- Satsnummer
 DA- Lotnummer
 FI- Eränumero
 EL- Αριθμός παρτίδας

EN- Recycle
 DE- Recyceln
 FR- Recycler
 NL- Recyclen
 IT- Riciclare
 ES- Recicla
 PT- Reciclar
 SV- Återvinna
 DA- Genbrug
 FI- Kierrätä
 EL- Ανακυκλώστε το

EN- Use by date
 DE- Verfallsdatum
 FR- Date de péremption
 NL- Uiterste gebruiksdatum
 IT- Utilizzare entro
 ES- Usar antes de
 PT- Usar até a data
 SV- Använd före datum
 DA- Udløbsdato
 FI- Viimeinen käyttöpäivä
 EL- Ημερομηνία λήξης

EN- Not for injection.
 DE- Nicht zur Injektion.
 FR- Ne pas injecter.
 NL- Niet injecteren.
 IT- Non per iniezione.
 ES- No debe inyectarse.
 PT- Não injetar.
 SV- Använd för injektion.
 DA- Ikke til injektion.
 FI- Ei injektokäyttöön.
 EL- Μην ενέψετε.

Keep out of reach of children.
 Shelf life 36 months.
 Store at room temperature between 5°C-27°C (41°F-80°F).

U.S. Federal law restricts this device to sale by or on the order of a state licensed dental professional.
 For immediate reorder and/or complete descriptions of Ultradent's product line, refer to Ultradent's catalog or call Toll Free 1-800-552-5512. Outside U.S. call (801) 572-4200 or visit our website at www.ultradent.com.

General Precautions:
 1. For Professional use only.
 2. Review instructions, precautions, and SDS before beginning treatment. Use only as directed.
 3. Keep products out of heat/sunlight.
 4. Avoid skin exposure to resins.
 5. Isolate strong chemicals to area of treatment.
 6. Confirm that patient has no known allergies to treatment materials.
 7. Test flow of materials from syringe and tip before using intraorally.
 8. Never force syringe plungers.
 9. Use Ultradent® Syringe Covers and/or clean and disinfect syringes between patients.

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