

# PermaFlo® Pink

**EN-Flowable Composite**  
**DE-Flüssiges Komposit**  
**FR-Ciment dentaire composite liquide**

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**FI** Isolate preparation with retraction cord and/or rubber dam. Etch 20 seconds.

Apply puddled coat of PQ1, air thin, light cure 20 seconds.

Adapt PermaFlo Pink to appropriate contour. Verify no overhanging margins are present.

Complete contouring with appropriate burs, cups, disks, etc. Finish with Jiffy HiShine.

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## PermaFlo® Pink Flowable Composite

### Description & Indications:

PermaFlo Pink is a light-cured, radiopaque, methacrylate-based flowable composite. PermaFlo Pink contains 1µm average particle size with narrow upper limit particle distribution. It is thixotropic, 68% filled.

### Indications:

Use PermaFlo Pink to mask gingival recession where root structure is exposed.

### Procedure:

1. Remove luer lock cap and attach a disposable Micro 20 gauge tip securely onto syringe.
2. Isolate preparation with cord and or rubber dam to control sulcular fluid.
3. Etch 20 seconds.
4. Apply dentin bonding agent (We recommend PQ1). Air thin, light cure 20 seconds.
5. Apply first layer of PermaFlo Pink 0.5 - 1 mm and light cure. Follow with 2 mm increments as needed. Light cure each resin layer 20 seconds, see Fig. 1.
6. Shape with appropriate finishing burs. Finish with strips, polishing cups, and wheels or disks as appropriate. We recommend Jiffy® Polishing system for shine, ending with Jiffy HiShine, see Fig. 2.

### Precautions:

1. Partially cured material will have deficient physical properties, including strength, wear resistance retention, marginal integrity and color stability.
2. Incompletely cured material could also cause pulpal irritation. Use only high-quality curing units check units periodically with a light meter.
3. The oxygen inhibition layer should not be removed or disturbed during incremental layering because it optimizes the bond to each succeeding layer.
4. Any resin can be sensitizing to skin or can cause allergic reaction if individual is already sensitized. Therefore, wash area of skin exposed to resin completely to remove material. If dermatitis, rash or other allergic reaction is observed, see a physician.
5. Syringe caps should be replaced following use to prevent polymerization (as composite resins are activated).


Fig. 1

Fig. 2